


Key Takeaways

- LIMS a good fit for your labs overall sample management
- DIMS streamlines order entry and results approval workflows
- Datagain guides customer via a LIS capabilities decision tree

Purpose of LIMS

- Increase Sample Volume Throughput
- Reduce Staff Time Managing LIS
- Automate Result Distribution
- Insure Accurately Reported Results



Lab Automation

1. Order Entry



2. Results Approval



3. History and Reporting



Brown County Adult Drug and DUI Court Lab					
1328 Ashburn Ave, W 212 1000 Salem, NC 28786 732.466.2000					
Client Test History Details					
Client ID	Accession #	Panel	Conducted by		
0254707-0000-0106-0002-FY-0000	00541310	75 Panel 7000	Move Smith		
Collection	Specimen	Agency	Northwestern Drug Court and 200-7000		
Client name	Blairford, Robert	Results Exposed	2017-03-02 15:02:36pm		
Results	Ref. range	Outcome	Collection date	Results date	Comments
Alcohol	< 0.10	Negative	2017-03-02 15:02:36pm	2017-03-02 15:02:36pm	Lab: 0000000000000000
Urine/serum	0.00-0.10	Negative	2017-03-02 15:02:36pm	2017-03-02 15:02:36pm	Lab: 0000000000000000
Urine/serum	0.00-0.10	Negative	2017-03-02 15:02:36pm	2017-03-02 15:02:36pm	Lab: 0000000000000000
Urine/serum	0.00-0.10	Negative	2017-03-02 15:02:36pm	2017-03-02 15:02:36pm	Lab: 0000000000000000




Data Accuracy

- Unique Client ID
- Facial Recognition - Optional
- Sample Barcode Generation
- Site Specific Fees & Payments



Imagine If

- Your lab was automated and efficient
- Integrated with Sure Digital Cup results
- Donor history is always accessible
- No mistakes made from check-in to results
- You were 100% paperless



Communicate & Collaborate More

- Specialty Court Case Management
- Color code randomization module
- Call-to-schedule module
- Probation & Pretrial Management App

Factors when considering a LIMS

- Customer Type – CJ, Tox, Ref, Treat
- Test Volume
- Lab Configuration – Instrument or Cup
- Speed to go live and ease of use
- Data sharing requirements
- Interface and Integration

Conclusion

- LIMS is an essential tool for labs
- DIMS LIS increases testing throughput in less time
- Being able to configure then control your lab is critical to success

Contact Us

 dims@datagainservices.com

 (606) 391 3467  <https://dimsservices.com>